

**BPL
TECHNOLOGY
FOR LIFE**



Believe in the Best



OXYGENATORS OG 4203

BPL OG 4203

- Shorter setup and testing time
- Maintenance free rotary poppet valve
- Dual Audio alarm system
- Portable, user-friendly, safe, quiet & reliable



FEATURES AND BENEFITS

- ABS injection moulded cabinet which out performs PVC and Noryl in strength and durability
- Flow Control Knob allows you the ability to protect your customers from the harmful effects of exceeding prescribed oxygen level
- 4.35 ~ 10.15 PSI oxygen output pressure allows the patient the freedom to use longer tubing without restricting the flow
- Quiet unit – allowing restful sleep for the patient
- No other raw materials are needed except air and low cost of oxygen production
- All gas circuit is controlled by low pressure system, which operates securely and reliably
- Adjustable flow rate; convenient operation; excellent performance
- Built-in atomising device for atomising therapy (Nebulisation therapy)
- Uninterrupted flow of oxygen

TECHNICAL SPECIFICATIONS

GENERAL	Flow Rate	Adjustable flow rate 0 ~ 5 LPM
	Purity	≥ 90%
	Oxygen delivery pressure	0.03 to 0.07 Mpa (4.35 ~ 10.15 PSI)
	Atomising pellet (ml / min)	≥ 0.5
	Alarms	Low Pressure High Pressure Power Failure
	Filters	Gross Particle Cabinet Filter Compressor Intake Filter Bacterial Filter
	Sound Level Noise	< 52 db (A)
	Electrical Requirements	220 ± 10 % VAC; 50 Hz... 2 amps
	Power Consumption	≤ 500 watts
	Size	640 mm (H) x 410 mm (W) x 410 mm (D)
Weight	32 kg.	
ENVIRONMENTAL	Operating temperature	0° C to 40° C
	Relative Humidity	≤ 85 %
STANDARD ACCESSORIES	Humidifier bottle	2 nos.
	Nasal Canula with extension tubing	2 nos.
	Operating manual	1 no.

* Due to constant upgradation, specifications & features subject to change without prior notice.

Marketing Office :

BPL Limited, Healthcare Business Group, BPL Towers, 13, Kasturba Road, Bangalore - 560 001. India. Ph: + 91-80-22270134. Fax: + 91-80-22275641

BPL Branch Sales Offices :

Ahmedabad	66314471	98240-51286
Bangalore	41144681	98866-45749
Bhopal	4262001	98260-20615
Cochin	3251718	98464-78287
Chennai	28113170 / 71	98400-55538
Chandigarh	4626572-74	98151-12908
Hyderabad	23261103 / 04	93924-38280
Jaipur	2681576	99283-68019
Kolkata	23648951 / 52	99030-11602
Lucknow	4042194 / 2439479	94154-18773
Mumbai	24173295 / 24173210	98700-00738
New Delhi	32915096 / 32915097	98101-00930
Patna		93349-01454

Your nearest BPL Channel Partner address



For enquiries contact: sales.healthcare@bplmail.com

www.bplworld.com/healthcare

