

## Babytherm® 8004/8010 Neonatal Open Care

The latest generation of Babytherm® infant warming systems from Dräger sets new standards in open care, not only in providing babies with a unique thermal environment, but also in ensuring you enjoy outstanding ease of access so you can give them the best possible care. Babytherm® is an intensive care Infant Warming System intended as a therapeutic device providing warmth for premature and full term neonates.



## Benefits

---

### Optimal Thermoregulation

Babytherm works to deliver an optimal thermo-neutral environment for the neonate safe from cold or heat stress. The Babytherm's ceramic heater elements and golden reflectors provide even heat distribution over the entire bed area. So wherever you place the baby on the bed, your little patient will always feel comfortably warm. The Babytherm also provides a fully integrated heated conductive gel mattress to keep the baby warm at all times. Measuring a baby's central and peripheral (toe) temperature gives you an early warning of thermal instability or possible cold stress by observing the measured values on the integrated digital display.

---

### Ease of Use

Our efforts to ensure you enjoy the best possible working conditions are reflected in a number of practical features. The control panel positioned directly beneath the radiant warmer is easy to reach so you can adjust the settings with no trouble at all. Another advantage is that you will easily be able to read all the essential information on the large, bright digital display from any position around the bed. The bright red central warning light at the end of the radiant warmer will alert you to an alarm condition even if you are working with another baby or in a different part of the room. A clear text message then provides you with all the details you require. The height of a Babytherm can be adjusted to suit your height, swivel cupboards are available to give you extra storage space for accessories and all the moving parts or control functions are easy to handle.

---

### System Synergy

Together with the intelligent design of architectural systems, bed side monitoring, and gentle neonatal ventilation, you enjoy an integrated work environment that is ergonomically smart, intuitively simple, and perfectly accessible for improved communication, time savings, and decreased risk of errors across the NICU and beyond.

---

### Intelligent Accessibility

The smartSWIVEL mechanism automatically focuses radiant warmth on the baby, even if the radiant heater is moved to one side to allow free working space for the caregiver. As a result, you don't need to worry about the baby cooling down or suffering cold stress during X-rays or other procedures because you can be sure of even radiant heat distribution at all times.

---

### Developmental Care

Babytherm's developmental care concept comes in the smartSWIVEL which allows closeness for the family without compromising warming of the infant. Long dual temperature sensor cables facilitate the practice of Kangaroo Care. The examination lamp features settings for day and night. Two brightness levels provide just the right lighting anytime. A central alarm light could be set up or a ramping alarm sound that is a low pitch melody set up away from the bed.

## Accessories



D-506837-2012

---

### Dräger VarioLux®

The Dräger VarioLux® offers cool, highly efficient, and variable intensity examination light specially designed for use in neonatal intensive care settings. Made for single-handed operation, it features adjustable color settings as well as a dimmer function.



D-7384-2009

---

### Neonatal Care Accessories

Dräger accessories for neonatal ventilation and thermotherapy ventilate gently and effectively, reduce stress and help promote the development of the newborn with a wide range of accessories designed specifically for use with the smallest of patients.

## Technical Data

BABYTHERM	8010	8004
<b>STANDARD FEATURES</b>		
SmartSWIVEL mechanism	•	•
Manual mode	•	•
Skin servo mode	•	•
Bright and soft light	•	•
Central Alarm	•	•
ThermoMonitoring™	•	•
Text messages	•	•
Bed tilt	•	•
Inner side walls (7 cm)	•	•
Heated gel mattress	•	
X-ray tray	•	•
<b>OPTIONS</b>		
Integrated phototherapy	•	•
Height adjustment	•	•
RS232 interface	•	•
Storage Drawers	•	•

Not all products, features, or services are for sale in all countries.  
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to [www.draeger.com/trademarks](http://www.draeger.com/trademarks) to find the current status.

**CORPORATE HEADQUARTERS**  
Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
[www.draeger.com](http://www.draeger.com)

**Manufacturer:**  
Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck, Germany

**REGION DACH**  
Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION EUROPE**  
Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION MIDDLE EAST, AFRICA**  
Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

**REGION ASIA PACIFIC**  
Dräger Singapore Pte. Ltd.  
25 International Business Park  
#04-20/21 German Centre  
Singapore 609916  
Tel +65 6308 9400  
Fax +65 6308 9401  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**REGION NORTH AMERICA**  
Dräger, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 800 4DRAGER  
(+1 800 437 2437)  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

**REGION CENTRAL AND SOUTH AMERICA**  
Dräger Panama Comercial  
S. de R.L.  
Complejo Business Park,  
V tower, 10th floor  
Panama City  
Tel +507 377 9100  
Fax +507 377 9130  
[contactcsa@draeger.com](mailto:contactcsa@draeger.com)

Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

