


**The Reliable Fixed X-Ray Imaging solution**



www.bplmedicaltechnologies.com



Happier Living Everyday



Follow us on

# Regulatory approvals are Pending. Not for Commercial Sale until Regulatory approvals are secured.




A BPL Medical Technologies Magazine

Powered by



Scan the product image above with BPL AR App to view the product video of BPL X Rad 100




Download BPL-AR App to see this catalogue come to life  
Download BPL Promise App to experience the latest in medical technology

**BPL XRAD<sub>300</sub> Technical Specifications**

Generator Output	24 kW as per ISI (IS: 7620-Part-I).
mA max	300 mA
mAs range	Range Upto 500 mAs
KVp range	RANGE 40 to 125 KVp in Steps of 1 KVp/Step
Radiography Technique	Three point Technique (KV, mA & mAs)
X-Ray Tube	TOSHIBA (rotating anode E-7239X)
Focal Spot	large 2 mm <sup>2</sup> , small 1 mm <sup>2</sup>
Collimator	Manually operated with light beam diaphragm ON/OFF From Control Panel
Control Panel	Digital Control Panel with Alphanumeric LCD (20x4)
Anatomical Programming (APR)	Yes
Other Features	Memory to retain last selected factor Pre-Set Factors Self Diagnostic Program
Stand	Fixed Stand (Floor to Ceiling)
Hand Switch	Dual Action
Power Supply	440 V (±10%) AC 50 Hz, Three Phase, 63 Amps, Line Resistance of < 0.4 Ohms
Compatible Tables	BPL Plano – Horizontal Table BPL Multio – Multi Position Table BPL Floato – Table with Floating Top
Optional	Vertical Bucky Stand

**A Global Medical Technology Company**



BPL X-Rad 300's Service is Powered by



Service Helpline  
**1800-425-2355**

**CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY**

BPL Medical Technologies Private Limited  
11th KM, Bannerghatta Road,  
Arakere, Bangalore - 560076  
CIN: U33110KA2012PTC067282  
Ph: +91 80 26484388/ 2648 4350  
Toll Free: 1800-4252355  
Website: www.bplmedicaltechnologies.com  
For Enquiries: sales.medical@bpl.in



Happier Living Everyday